GLOBUS

LASER THERAPY

 *User Manual*

 **PHYSIOLASER 1000**

 **PHYSIOLASER 500**

CE
0476

DEAR CUSTOMER

THANK YOU FOR CHOOSING A GLOBUS PRODUCT. WE REMAIN AT YOUR ENTIRE DISPOSAL FOR ANY ASSISTENCE OR ADVICE YOU MAY NEED

Pysiolaser 1000

Pysiolaser 500

The lasertherapy devices (**Pysiolaser 1000 e Pysiolaser 500**) are produced and distributed by:

DOMINO s.r.l.

via Vittorio Veneto 52

31013 - Codognè - TV - Italy

Phone 0438.7933

Fax. 0438793363 www.globusitalia.com

info@globusitalia.com

These products have been manufactured according to the technical regulations in force and are certified according to Directive 93/42/EEC updated by 2007/47 directive for medical devices, by Kiwa CERMET s.p.a. Body (authorization n. 0476), in order to ensure the products' safety.

Table of contents

TECHNICAL FEATURES	6
Device	6
Use conditions.....	6
Technical features of the laser.....	6
Accessories	7
EQUIPMENT	8
CONNECTIONS.....	9
Device	9
Handpiece Laser	9
Battery: how to charge the batteries.....	9
LABELLING AND SYMBOLS.....	10
Device	11
Panel and keyboard.....	12
Display/Interface	12
ALARMS	13
Compliance	13
Meaning	13
Warning sounds	13
WARNINGS AND CONTRAINDICATIONS.....	14
Mandatory behavior.....	14
Warnings before the use.....	14
Warnings during the use	15
Side effects and contraindications	15
Contraindications.....	16
MAINTENANCE AND CLEANING.....	16
Maintenance and cleaning of the device	16
Battery	17
Disposal of the device	17
INSTRUCTIONS FOR USE	18
Pin unlocking.....	18
Program selection from the "Program List" menu	18
Visualization during a program execution	20
How to change the working phase parameters	21
Special Functions.....	21
Treatments.....	21
Activation of the STIM LOCK function	21
The "Programming" function.....	22
Scanning mode	22
Multi spot mode.....	23

Advanced Menu.....	24
APPLICATION MODE: THERAPEUTIC AND BEAUTY EFFECTS	26
WARRANTY CONDITIONS	30

TECHNICAL FEATURES

Device

Size: 160x99x35,4

Weight: 404 g.

Case: in food grade ABS

Protection level of the case: IP 22

Storage and transportation temperature: from -10°C to 45°C

Max. relative humidity: from 30% to 75%

Use conditions

Temperature: from 0°C to 35°C

Max. relative humidity: 15% to 93%

Atmospheric pressure: from 700 hPa to 1060 hPa

Technical features of the laser

Medical Device of II b class

Laser of IV Class

Nominal power (CEI EN 60825-1): Up to 1000mW

Beam divergence: 13°

Power Density: up to 1 W/cm² max (Physiolaser 1000)
500mW/cm² max (Physiolaser 500)

NOHD: 0.27 m (1W)

0.18 m (500mW)

Spot size: 3.5 mm or 11.5 mm

Wave length: from 630 to 1064 nm 808 nm

(the wave length of the laser source is indicated on the handpiece)

Frequency: 8 – 10.000Hz

Operation mode: Continuous and Pulsed

Wave length of the guidance light (I class device): 660 nm

Accessories

Glasses for the patient (Physiolaser 1000)

Visual field PPE: >40

Length regulation and sidepiece inclination: jerky mechanism

Ventilation system: small holes in the lateral protections

Complete PPE weight: 34

Materials: polycarbonate

VLT%: 0,36

Color: green 7

Bend: Basic 6

Optical class: 2

Scale number: 7

Glasses for the operator

Visual field: >40

Frame weight: 28 g

Frame material: polyamide 6.6 varnished

Transmittance (%): 35

Bend: Basic 6

Optical class: 1

Minimum thickness: 2.0

Power supply unit

Brand GLOBTEK

model: GTM41060

PRI: 100-240V~ 50-60 Hz 600 mA

SEC: 12 V = 2,08 A


Battery

Battery pack:

Ni-MH 7,2 V 1,8 Ah

EQUIPMENT

- A:** CARRYING CASE
- B:** GLASSES FOR THE OPERATOR
- C:** GLASSES FOR THE PATIENT
- D:** WARRANTY
- E:** POWER SUPPLY UNIT
- F:** LASER HANDPIECE
- G:** DEVICE

The device is supplied complete of laser handpiece, power supply unit, glasses for the operator and glasses for the patient; opening the package, it is necessary to check that the basic equipment is complete. If some elements should be missing, contact immediately the authorized retailer where you purchased the product. Control carefully the integrity of the device and its accessories.

INTENDED USE

The after sales service is guaranteed for 5 years. We suggest having a check of the device every two years for the maintenance and to ensure the safety.

The Physiolaser is designed to be used in the following operating environments:

- home environment;
- clinics;
- physiotherapy centers;

- rehabilitation centers;
- general pain treatments (medical treatments);
- for beauty and sport purposes.

CONNECTIONS

Device

Power supply by electricity grid

The Physiolaser operates connected to the mains. To connect the power supply unit to the connector, plug it as shown below.

Handpiece Laser

To connect the cable of the handpiece to the device plug it into the appropriate inlet removing the sticker with the yellow symbol (see picture).

Battery: how to charge the batteries


The device is supplied with a set of rechargeable nickel-metal hydrate batteries (7.2V, 1.5Ah), which have high performance without storage effect.

Recharge the batteries when the battery light on the display indicates $\frac{1}{4}$. To charge the batteries, turn off the device, then connect it to the power supply unit included by inserting the plug in the appropriate inlet (see picture).


Do not use a power supply unit different from the one provided with the device. To replace the device's batteries, contact an authorized service center.

LABELLING AND SYMBOLS

	Warning
	This symbol on your device indicates that it complies with the directives on medical devices (93/42/CEE 47/2007CEE). The number of the notified unit is 0476.
	It indicates that this is a I class device.
	It indicates that this device has BF type parts.
	<p>WEEE symbol (Waste of Electrical and Electronic Equipment). Recycling symbol.</p> <p>The WEEE symbol used for this product indicates that the device may not be disposed of as a household product. Properly dispose of the product to help protect the environment. For more information on recycling this product, contact the local competent department, the household waste management company or the store in which the product was purchased.</p>
	It indicates that the product has been produced respecting the directive 2002/95/CE.
	It indicates the optimal temperatures for the storage and transportation of the product.
	It informs the operator that before using the device he must read the manual.
IP22	It indicates the water protection degree
	It informs the operator of a compulsory conduct
	It refers to the pressure of the environment where the device and its accessories are stored and shipped
	It refers to the humidity of the environment where the device and its accessories are stored and shipped
	It refers to the manufacturer
	It refers to the expiry date


	It refers to the production lot
RH	It refers to the percentage of storage humidity
	It refers to the manufacturing date
	It refers to the interruption of the laser emission
	Attention danger laser emission not visible
	It identifies the risks related to the laser use and it classifies them
PRI	Mains voltage
SEC	Power supply voltage
Input	Input: It indicates the value of the mains voltage for the power supply unit
Output	Output: - it indicates the power supply unit outbound voltage - it indicates the maximum power value of the magnetic field emitted by the device - it indicates the range of frequencies of the magnetic field emitted by the device
Type	It indicates the device type
Power	It indicates the power supply unit model of the device
Battery	It indicates the battery pack inside the device

Device

Input	GTM4160-2512
PRI:	100-240Vac 50-60Hz Max 0,6A
Sec:	12VDC 2.08A
	12Vd.c. ; 0,55A 6,6VA
Output	808nm 1000mW
Laser Type:	LAS
LAS051309199	
	LAS051309199
	Batteries Ni-Mh 7.2V 1800mAh

Panel and keyboard

ON/OFF/OK Button (To turn on/off and confirm the treatment)

LEFT/BACK Button

P+ Button (to scroll the programs up)

P- Button (to scroll the programs down)

RIGHT Button

Up/Down left button to modify the power

Fn button to activate some special functions described in the "instructions for use"

Display/Interface

Sounds disabled

Program name	Set Prog. Libero	Battery
Spot diameter	Ø: 3.5mm	Handpiece animated icon
Wave length	λ: 0nm	NO CABLE
Frequency	F: 1000Hz	Energy per point
Duty cycle	d: 50%	Points treated/ to be treated
Treatment area	A: 50cm ²	Emitted energy
Dosage	D: 0.1J/cm ²	Peak power
Remaining time of the program	T: 1min 0sec	

Additional display elements: *




ALARMS

Compliance

Certifications: CE MDD Certificate

The sound and acoustic signals are in compliance with the directive 60601-1-8.

Meaning

	<u>Correct power emission:</u> the treatment proceeds correctly, the laser emission is present.
	<u>Cable not connected</u> Verify that the cable is connected to the right output.
	<u>Device in pause.</u> Press OK to continue the treatment.
	<u>The device is in stand-by mode: Press the button on the handpiece to start the emission</u>

Warning sounds

If the Warning Sounds function is on, the device emits an acoustic signal while pressing the buttons.

If the Warning Sounds function is off, the acoustic signal that indicates the start and the end of the treatment and the start and the stop of the laser emission are still present.

If during the treatment the laser handpiece is suddenly disconnected from the connector, the device will emit an acoustic signal.

WARNINGS AND CONTRAINDICATIONS

Mandatory behavior

For safety reasons, the device must be used in the prescribed manner and within the limits of use explained in this manual.

The manufacturer declines any responsibility due to a different use from what it is indicated in this manual.

No part of this manual (texts and/or photos) may be reproduced by electronic or mechanical means without the prior written authorization of the manufacturer.

If the package, the cable or the connector of the power supply unit show signs of wear or damage, replace them instantly.

The unit should be connected to the mains by its power supply unit. Before starting the treatment, make sure that the power system specifications comply with the directives in force within your country. Make sure that the power supply unit will be in a comfortable position and that it will be easy to be removed.

The device should be protected against a not authorized use.

It is forbidden to modify the device.

Warnings before the use

We suggest avoiding the use of the device together with other electronic devices, especially with those used to maintain vital functions. For a correct use of the Electromedical device, make reference to the tables attached. If it is necessary to use the device near or together with other devices, pay attention to its working.

- It is recommended to read carefully the entire operating manual before using the device; keep carefully this manual.

- Do not insert fingers or metallic objects into the emission holes to avoid injuries and/or damages to the optical section.

- The unit should be used only with the laser handpiece provided originally (or Globus branded) and according to the therapeutic modalities described.

- Before each use always check the integrity of the unit. This is a fundamental requirement for carrying out therapies, do not use the unit if either the buttons or cables are defective or malfunctioning.

- It is recommended to control the functioning of the pointing laser beam. This must always be present before executing the treatment; if it should be missing, it means that an anomaly occurs.

The device:

- should not be used in front of people who are not wearing the protection glasses.

- it must be used only from people over 18 who are able to understand and take action.

- it must be used following the indications and under the physician or physiotherapist's control.
 - it must be kept out of the reach of children.
 - if there is any health problem, it must be used only after seeing a doctor;
 - once you have turned on the device, make sure that in the display the software version and the model of the device appear, it means that the device is working and it is ready to be used.
- If not, or in the display all the segments do not appear, turn it off and on again. If the problem persists, contact the service center and do not use the device.
- An unexpected switching-off of the device means that the battery has run-down. Charge it according to the instructions in the paragraph: HOW TO CHARGE THE BATTERIES
 - The use of controls or the execution of different procedures from what specified can cause a dangerous exposition to the laser radiation.

Warnings during the use

During the use of the Physiolaser some warnings should be followed:

- the device should be used wearing the included appropriate glasses for the operator and for the patient.
- In case of damaged cables, they must be replaced with original parts and not used anymore.
- The device must be kept out of the reach of any pet that could damage it and contaminate it with parasites.
- The cables of the laser handpiece and of the power supply unit should not be wrapped around people's neck to avoid any risk of strangulation and suffocation.
- The mobile and fixed radiocommunications devices could influence the functioning of the electromedical device: see the tables attached to this manual.
- The use of controls or the execution of different procedures from what specified can cause a dangerous exposition to the laser radiation.

Side effects and contraindications

Possible side effects of the lasertherapy are:

- photosensitization;
- rebound effect (the increase of the tenderness after the lasertherapy use);
- skin pigmentation.

Contraindications

Do not use the lasertherapy in the following cases:

- Do not carry out this treatment on pregnant women.
- On children under 15.
- Do not point the laser on big moles or on dark areas that could absorb excessive quantity of energy. In this case we suggest covering these areas with a titanium dioxide or zinc cream.
- On noble organs (ovaries, heart, head, etc.).
- On epileptic patients.
- In case of fibrocystic mastitis.
- On patients who have photosensitization reaction.
- In case of fever.
- In case of patients with tumor diseases, see your physician or oncologist before using the lasertherapy.
- Do not irradiate the thyroid and other endocrine glands.

MAINTENANCE AND CLEANING

Maintenance and cleaning of the device

- In case of real or alleged malfunction, do not tamper with the device or try to repair it by yourself.

Do not intervene on the device and do not open it. Only specialized and authorized centers can repair it.

- Avoid violent impacts that may cause damage and malfunctions to the device even if undetectable immediately.
- Use this device in a dry environment and in an open space (not wrapped in any materials).
- Clean the device and its accessories only with disinfectant with sodium hypochlorite or quaternary ammonium salt diluted with distilled water equal to 0,2-0,3%. After cleaning/disinfecting it, dry the device and its accessories with a clean cloth.
- Use always the device and its accessories with clean hands.
- It is recommended to use the device in a clean room, to avoid the contamination of the device with dust and dirt.
- It is recommended to use the device in a ventilated space, with regular air change.
- It is recommended to clean/disinfect the parts after every use.

Battery

Battery info.

The device is equipped with a menu that can visualize the charge of the battery, the values and the conditions of the end of the charge.

It is recommended to access to this menu only after a complete charge of the batteries.

From the main menu choose "Advanced", then "Set-Up" and finally "Battery info".

Six codes will be visualized with the following meanings:

COD1 = 0 expected voltage threshold reached.

COD1 = 1 max. charge time reached.

COD2 value of the battery voltage at the start of the charge.

COD3 value of the battery voltage at the end of the charge.

COD4 charge duration (from 1 to 840 minutes, ideal time 720 minutes).

COD5 charger/power supply unit connection duration.

COD6 Battery pack voltage value.

According to the values above described, it is recommended to replace the battery when COD1 = 1 and COD3 < 7,4 volts. Or when COD3-COD2 \geq 2 volts and COD4 < 600. Or when COD6 is lower than 5,8 volts.

Furthermore it is recommended to replace the battery pack after 3 months in which the device has not been used. After that period, batteries usually loose their ability to charge making the recharge dangerous.

Disposal of the device

Do not throw the device or part of it in the fire, but dispose of the product in the specialized centers and respecting the directives in force within your country. When the product has to be disposed of, the user can give it back to the retailer when purchasing a new unit.

A correct separate waste collection or following what above mentioned contributes to avoid possible negative effects on environment and health and promotes the reuse and/or recycle of materials of which the device is composed. The illegal disposal of the product by the user involves the application of the administrative fines according to the current regulations.

INSTRUCTIONS FOR USE

Pin unlocking

The device, when not used, should be protected against a not authorized use. The PIN function has this objective; it consists of a SW code intended to impede the device use to not authorized staff.

For a correct connection of the product it is necessary to proceed as follow:

- Connect the handpiece to the relative output of the device;
- To turn on the device, press and hold the ON/OFF (OK) button for approximately 3 seconds until a sound signal is heard.

Pin unlocking procedure

Press 4 times the P+ or "Save" buttons.

To switch-off, press and hold the ON/OFF (OK) button until a sound signal is heard.

With the Up and Down buttons of the joypad, you can move inside the main Menu, on the following entries: "Last 10", "Program List", "Treatments", "Programing", "Advanced", "Switch-off".

By pressing the OK button on "Last 10" you access to the Menu: the last 10 executed programs are shown.

Program selection from the "Program List" menu

Use the UP and DOWN buttons to place the cursor on "Program List" and confirm by pressing the OK button.

Area selection

By using the Up and Down buttons, place the cursor on the area you want to select. Press OK to confirm.

Program selection

Select the desired program by using the Up and Down buttons.

Press OK to confirm.

Press and hold the left (Back) button for 3 seconds to return to the previous screen.

Select START to start the treatment.

Select SAVE IN TREATMENTS to save the program and execute it later (see paragraph about "Treatments").

Phototype selection

Once you have chosen the program, you were asked to set the skin phototype of the subject to treat:

Phototype setting
Light (A)
Dark (B)


By choosing the phototype B, the emission intensity will be automatically reduced and the duration of the treatment will be increased, both of 20 %.

Light phototype= for light/normal skin


Dark phototype= for dark skins or dark marks.

How to start the program

After selecting the program, the display will show the message "ready". At this point put the handpiece on the area to be treated and press and hold the emission button . The handpiece emits only if the button has been pressed, if not, the device does not emit.


Visualization during a program execution

Visualization of the parameters during the execution

By holding the * button, the screen will show you the visualization normally characterizing a lasertherapy protocol. The following screen with the specific values of the ongoing program will appear:

How to pause/stop the program

To pause the program, press the OK button. Press OK again to return to the program. During the pause the message PAUSE will appear on the screen.

How to exit from the program

To stop the program before the end of the treatment, press the Back button.

How to change the working phase parameters

After starting the program, the power of the treatment can be modified by using the P+ and P- buttons. By increasing or decreasing the power, the duration of the treatment will change.

Special Functions

“Last 10” Function

The device stores the last 10 executed programs, so that these are available for a rapid and easy execution.

The storage occurs automatically at the end of each program; when the memory is full, older programs are automatically deleted.

When the device turns on, select "Last 10" and then confirm with OK

Treatments

The device allows the user to block the device and to execute only the treatments saved with the “Save in treatments” function in the screen previously to the execution of the same program. This function is especially thought for the rent of the unit to inexperienced users and/or to patients that have to carry out only some special programs that have been chosen by the specialist.

Activation of the STIM LOCK function

Press and hold the buttons fn and --> (RIGHT button) for at least 3 seconds until the area where the treatments have been saved appears. After the activation of the STIM LOCK function, the functions of the unit will have limited functionalities.

Deactivation of the STIM LOCK function

Press and hold the fn and <-- buttons (LEFT button) for at least 3 seconds until the main menu appears.

NOTE: If, when the unit is turned on, the main menu does not appear, verify that the STIM LOCK function is not activated. Try to deactivate it.

If the problem persists, contact the customer service.

The "Programming" function

The "Programming" function is activated only when the laser pointer is connected. The device offers the possibility to create/modify new programs. This makes the device flexible and adaptable to all users' requirements.

From the "Programming" entry of the main menu, it is possible to create new programs (when the message EMPTY appears) and to execute already personalized programs. These programs can be modified at any time. The programs created in the "Programming" menu are the same for all USERS and cannot be stored in the "Last 10" list.

How to create a new program

From the main menu, move the cursor on "Programming" using the UP and DOWN buttons and confirm with OK.

Use the DOWN button to scroll down the memory list (from 1 to 15) and confirm with OK. The free memories are the ones where the message "EMPTY" appears.

NOTE: If the program has already been created, see "How to modify a program".

Insertion of the program's name

Use the LEFT and RIGHT buttons to select letters and confirm with OK. To delete a letter, move the cursor on DELETE and confirm with OK.

Parameter insertion


After selecting the program name you should select the desired carrier:

- Continuous (duty cycle 100%)
- Pulsed (duty cycle 50%)

After that, you were asked to choose to execute a program in multi spot mode or in scanning mode.

Scanning mode


When you choose the scanning mode, the following screen appears:

By using the Up and Down buttons you can modify the frequency (only in pulsed mode), the program duration, the intensity and the size of the area to be treated, and you can combine these parameters in order to have the desired protocol. In the upper part of the screen both the emitted total dosage and the dosage per cm² will appear.

Multi spot mode

When you choose the multi spot mode, the following screen appears:

By using the regulation buttons you should set the desired intensity for the treatment, by choosing the point number and the dosage per point.

By pressing the * button you can visualize and modify the information "E/p" Energy per point or the information "D (dosage) in J/cm²".

This allows the operator to choose the parameter he wants to set.

How to modify a program

Inside the "Programming" menu, select the program number you wish to modify and press at the same time the "fn" + UP" buttons.

How to delete a program

Inside the "Programming" menu, select the created program you wish to delete and press at the same time the "fn" + Down" buttons.

The deleted program cannot be recovered.

Advanced Menu

By selecting the Advanced Menu you can access to an area intended for the special functions.

"User Selection"

The "User Selection" in the advanced menu function permits to use the special functions ("Last 10 " and "Treatments") in a personalized manner.

NOTE: Every time the device is turned on, the last user will be displayed.

SETUP

- Backlight time FUNCTION

The backlight time function in the setup menu allows to modify the duration of the light during the stand-by phases using the Up and Down buttons.

- Contrast FUNCTION

The contrast function in the setup menu allows to modify the display contrast by pressing the UP/DOWN buttons.

- Auto shutdown FUNCTION

The "Auto shutdown" function in the setup menu permits to choose, using the UP and DOWN buttons, a period of time (in minutes) after which the device, if not used, automatically turns off.

- **Language selection FUNCTION**

The "Language Selection" function in the setup menu permits to choose among 4 different languages for the navigation using the UP and DOWN buttons. Confirm the selection with OK.

- **Warning sounds FUNCTION**

The "Warning sounds" function in the setup menu permits to enable (ON) or disable (OFF) the acoustic tones emitted by the unit.

- **"Handpiece button" FUNCTION**

The handpiece button function permits to select two modes (monostable and bistable).

Monostable: The device emits while pressing and holding the handpiece button; the device does not emit releasing the button.

Bistable: By pressing once the handpiece button the device emits, and by pressing another time the button the device does not emit.

- **"PIN change" FUNCTION**

The "PIN change" function present in the Setup area permits to modify the current PIN.

- Insert the current pin and press OK.
- Insert the new pin and press OK.
- Confirm the new pin and press OK.

The new pin is enabled.

- **"Diameter spot laser" FUNCTION**

In this section it is possible to set the diameter spot of the handpiece to use during the treatment. This value is indicated in the table with the program list in the end of the manual.

When choosing the 3.5 diameter you should remove the lens from the handpiece unscrewing it in a counterclockwise sense.

- **"Battery info" FUNCTION (see the chapter "MAINTENANCE AND CLEANING")**

APPLICATION MODE: THERAPEUTIC AND BEAUTY EFFECTS

The biological changes caused by the laser light, determine the following effects:

- 1) *Biostimulating effect.* The biostimulating effect of the laser accelerates the cicatrization of ulcers and sores. The laser light stimulates the mitochondris and accelerates the ATP production. The increased production of APT recharges the cells of energy so that, if the cells are damaged because of inflammatory, traumatic or regenerative causes, they restart to execute their physiological functions.
- 2) *Activation of microcirculation.* The lasertherapy has an intense vasoactive action on the microcirculation. The activation of microcirculation favors a greater nutritional contribution and a better drainage of catabolites from the tissues.
- 3) *Antalgic effect.* The laser has an analgesic effect because it increases the excitement of receptors and it creates an antiphlogistic action.

Application method

To carry out lasertherapy sessions it is necessary to follow these indications:

- 1) It is obligatory for the operator and the patient to wear protective glasses for the retina and to install the device in environments without reflecting surfaces.
- 2) The body area should not be covered with clothes and the skin has to be cleaned with alcohol. The presence of superficial fat could cause reflecting phenomena of beam obstructing its penetration.
- 3) During the treatment the handpiece must be put in contact with the skin and the ray must have a 90° of incidence. The more you go away from the skin, the bigger the treatment area is and the less the power density is. Only the operator can use the handpiece unit. The laser handpiece must be moved with a circular movement or with a grid scansion movement.
- 4) The treatment time for each cm² of surface changes in relation to the type of used laser, to the emission methods (continuous or pulsed) and to the average power.
- 5) The total number of therapeutic sessions is correlated to the type of pathology to be treated (acute or chronic); it can vary from 10 to 20 sessions.

MAIN AREAS TO BE TREATED

SHOULDER PAIN		ANKLE SPRAIN
		
KNEE ARTHROSIS/OSTEOARTHRITIS		FINGER PAIN
		
ACHILLES TENDON INFLAM.		EPICONDYLITIS
		
BACK PAIN WITH RADICULOPATHIES		NECK PAIN
		
CARPAL TUNNEL		TEMPORO/MANDIBULAR PAIN
		

Program list HEALTH- REHABILITATION

In compliance with the last regulations of the Ministry of Health, medical treatments are indicated in the below list

The following programs are medical

Health	Physiolaser 500	Physiolaser 1000	Session/ week	# weeks	Ø Spot mm
Arthrosis	x	x	5	2	11,5
Ankle sprain	x	x	5	2	11,5
Knee pain	x	x	5	2	11,5
Ankle sprain edema	x	x	5	2	11,5
Meniscus pain		x	5	2	11,5
Knee osteoarthritis		x	2	4	11,5
Achilles tend inflam.		x	5	3	11,5
Finger pain		x	3	3	11,5
Myofacial shoulder pain		x	5	2	11,5
Articular shoulder pain	x	x	5	2	11,5
Epicondylitis	x	x	2	4	11,5
Frozen shoulder		x	2	4	11,5
Carpal tunnel		x	5	2	11,5
Chronic back pain	x	x	5	2	11,5
Acute back pain		x	5	2	11,5
Back pain with radiculopathies		x	5	4	11,5
Neck pain	x	x	2	7	11,5
Neck pain with radiculopathies		x	5	3	11,5
Myofacial neck pain		x	5	2	11,5
Neck osteoarthritis		x	5	2	11,5
Temporo- Mandibular pain		x	2	5	11,5
Tendinopathies	x	x	5	2	11,5
Injuries		x	2	5	11,5

Because of the presence of clinical programs, this product is a medical device. Therefore it is certified by the Kiwa CERMET s.p.a n°. 0476 according to 93/42/CEE directive for medical devices. The certification covers clinical applications.

Program list BEAUTY – FITNESS (non medical treatments)

X = Available program

Beauty	Physiolaser 500	Physiolaser 1000
FACE BEAUTY		
Wrinkles	x	x
Acne	x	x
Dehydrated skin		x
Skin rejuvenation	x	x
Skin – face lifting	x	x
Expression wrinkles	x	x
Neck wrinkles		x
Couperosis	x	x
Telangiectasias		x
BODY BEAUTY		
Breast stretch marks	x	x
Thighs, hips and gluteus edematous cellulite	x	x
Thighs, hips and gluteus fibrous cellulite		x
Body fat and spread cellulite	x	x
Stretch marks	x	x
Epilation		
Scares	x	x
Loc. lymphatic stasis		x

CE0476 does not refer to non-medical treatments.

WARRANTY CONDITIONS

The device is guaranteed to the first user for a period of twenty-four (24) months from the date of purchase against defects in materials or of the manufacturing, provided that it is used properly and maintained under normal operating conditions.

In order to take advantage of the warranty service, the user must comply with the following warranty clauses:

1. The products, and all accessories, to be repaired must be sent by and at the expenses of the customer in their original packages.
2. The product's warranty is subject to the production of a fiscal document (fiscal receipt, receipted bill or sale invoice), attesting the product's purchase date.
3. The repair work shall have no effect on the original warranty expiry date and shall neither renew nor extend it.
4. If no defects are found, during the repair work, the costs related to inspection times shall in any case be charged.
5. The warranty becomes void where the fault has been caused by: impacts, falls, erroneous or improper use of the product, use of non-original power supplies or charger, accidental events, alteration, replacement/detachment of the warranty seals and/or tampering with the product. The warranty does not cover damages caused during transportation when unsuitable packages are used (see point 1).
6. The warranty does not cover the inability to use the product, other incidental or consequent costs or other expenses incurred by the purchaser.

NOTE: Before returning the device for repairs, we recommend to read carefully the user instructions contained in the manual.

In case of any faults caused to the device, refer to your retailer.

The manufacturer reserves the right to make all necessary modifications at any time in order to improve the aspect and quality of the product.

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GUIDANCE AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC IMMUNITY – FOR ALL EQUIPMENT AND SYSTEMS

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The PHYSIOLASER 1000 is intended for use in the electromagnetic environment specified below. The customer or the user of the PHYSIOLASER 1000 should assure that it is used in such an environment.

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GUIDANCE AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC IMMUNITY – FOR EQUIPMENT AND SYSTEMS THAT ARE NOT LIFE-SUPPORTING

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RECOMMENDED SEPARATION DISTANCES BETWEEN PORTABLE AND MOBILE RF COMMUNICATIONS EQUIPMENT AND THE PHYSIOLASER 1000 FOR EQUIPMENT AND SYSTEM THAT ARE NOT LIFE-SUPPORTING

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The PHYSIOLASER 1000 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the PHYSIOLASER 1000 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the PHYSIOLASER 1000 as recommended below, according to the maximum output power of the communications equipment.

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For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

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